



1000 Western Drive
Brunswick, OH 44212
PHONE: 330-225-3200
FAX: 330-225-1499
www.columbiachemical.com

COLCLEAN SOAK 70

TECHNICAL DATA
05-03-2021

COLCLEAN SOAK 70 STEEL SOAK CLEANER

| | |
|------------------|--|
| COLCLEAN SOAK 70 | is a non-chelated, free-flowing granular blend of cleaning agents, builders, and surfactants. |
| COLCLEAN SOAK 70 | specially developed multi-metal cleaner to handle softer metals like copper, brass, and zinc die-cast. |
| COLCLEAN SOAK 70 | may be used on steel substrates and contains phosphates. |
| COLCLEAN SOAK 70 | incorporates a unique wetting system, providing excellent detergency and emulsifying properties. |

OPERATING PARAMETERS

| | |
|------------------|----------------------|
| COLCLEAN SOAK 70 | 8 - 10 oz per gallon |
| Temperature | 170 - 200° F |

ANALYTICAL PROCEDURE

1. Pipette 10.0 mL sample into a 250 mL Erlenmeyer flask.
2. Add 3 drops phenolphthalein.
3. Titrate with 1.0N HCl to color change from pink to clear.

Concentration (opg or % by volume) = mL Titrant • 2.25

HANDLING & STORAGE

Columbia Chemical recommends referring to the specific product Safety Data Sheets for safety, handling, and storage precautions.

NON-WARRANTY

The data contained in this bulletin is believed by Columbia Chemical Corp. to be accurate, true, and complete. Since, however, final methods of use of this product are in the hands of the customer and beyond our control, we cannot guarantee that the customer will obtain the results described in this bulletin, nor can we assume responsibility of the use of this product by the customer in any process which may infringe the patents of third parties.

